

Key Features

- Lightweight design allowing for easy patient application and comfort
- Fully encapsulated components for easy routine cleansing
- Specially designed to minimise artefact caused by patient movement, ambient light and electrical interference
- Cables use high quality alloy conductors with KEVLAR reinforcement and 100% medical grade silicone jacket
- Latex Free



SoftTip probe large



SoftTip probe small



SoftTip probe medium

Device Compatibility

Embla/ Embletta (Nonin oximeters with 9way D type connectors)

- Large 8000SL (max finger diameter app 25 mm)
- Medium 8000SM (max finger diameter app 20 mm)

Konica Minolta for Pulsox 3/ 3i/ 3iA series and Pulsox 300/300i series

- Large Env-L-KM (max finger diameter app 25 mm)
- Medium Env-M-KM (max finger diameter app 20 mm)
- Small Env-S-KM (max finger diameter app 15 mm)

These sensors are not available for Masimo oximeters (used in Visi/ Grey Flash)

Contact Information

For more information regarding the SoftTip® Sensors and compatibility and use with Stowood products, please contact us per details provided below:

Stowood Scientific Instruments

Beckley, Oxford, OX3 9UP
Email: sales@stowood.co.uk
Web: www.stowood.com
Tel/fax: +44 (0) 18653 58860

Also, refer to www.envitec.com for more information on Envitec, Denmark.